

absorbable

Poly-Dox™

Violet, monofilament, synthetic absorbable suture composed of polydioxanone

comparable to PDS II™

Benefits

- Extended holding strength
- Excellent handling with exceptional suppleness
- Good knot security
- Non-wicking helps reduce bacterial migration along suture line
- FDA-approved D&C Violet No. 2 color for better visibility during surgery
- Also available in Antibacterial PLUS line

Uses

- Extended wound support
- Soft tissue repair including muscle, stomach, intestine, uterus, bladder, ligature
- General wound closure

Wound Support

75% at 14 days
66% at 28 days
53% at 42 days

Absorption

180 days



The RIGHT products at the RIGHT price.



absorbable

Poly-Dox™

Violet, monofilament, synthetic absorbable suture composed of polydioxanone

comparable to PDS II™



POLY-DOX

Needle Type	Packaging Style	Code	Size	Needle	Length	Ethicon PDS
PHOENIX SELECT						
▽	⊖	PS0467	0	PCP-1	30"	Z467H
▽	⊖	PS0970	0	PCP-2	30"	Z970H
▽	⊖	PS0468	1	PCP-1	30"	Z468H
▽	⊖	PS0466	2/0	PCP-1	30"	Z466H
▽	⊖	PS0969	2/0	PCP-2	30"	Z969H
▽	⊖	PS0451	2/0	PFS-1	30"	Z451H
▽	⊖	PS0968	3/0	PCP-2	30"	Z968H
▽	⊖	PS0452	3/0	PFS-1	30"	Z452H
▽	⊖	PS0398	3/0	PFS-2	30"	Z398H
▽	⊖	PS0397	4/0	PFS-2	30"	Z397H
⊙	⊖	PS0340	0	PCT-1	30"	Z340H
⊙	⊖	PS0334	0	PCT-2	30"	Z334H
⊙	⊖	PS0317	2/0	PSH	30"	Z317H
⊙	⊖	PS0339	2/0	PCT-1	30"	Z339H
⊙	⊖	PS0333	2/0	PCT-2	30"	Z333H
⊙	⊖	PS0316	3/0	PSH	30"	Z316H
⊙	⊖	PS0332	3/0	PCT-2	30"	Z332H
⊙	⊖	PS0315	4/0	PSH	30"	Z315H

POLY-DOX

Needle Type	Packaging Style	Code	Size	Needle	Length	Ethicon PDS
⊙	∞	0338	3/0	NCT-1	30"	Z338H
⊙	∞	0311	3/0	NSH-1	30"	Z311H
ANTIBACTERIAL						
▽	∞	PM451	2/0	NFS-1	30"	PDP443H
▽	∞	PM452	3/0	NFS-1	30"	PDP442H
⊙	∞	PM316	3/0	NSH	30"	PDP311H
CASSETTES						
N/A	Cassette	SMC10	#0	N/A	25M	
N/A	Cassette	SMC01	#1	N/A	25M	
N/A	Cassette	SMC02	#2	N/A	25M	
N/A	Cassette	SMC20	2/0	N/A	25M	
N/A	Cassette	SMC30	3/0	N/A	25M	
N/A	Cassette	SMC40	4/0	N/A	25M	

⊙ Taper Point ▽ Reverse Cutting ▽ Precision Reverse Cutting ∞ Suture is wound in a Figure 8 package ⊖ Suture is wound in a Racetrack package



Needle silhouettes are for display purposes only